Sepax Hydrophilic Interaction Chromatography (HILIC) Separation – More Choices for Your Compounds

Innovative Surface Chemistry

Column Information

To solve the challenges of more and more highly polar pharmaceuticals and small biological molecules, Sepax developed a series of chemistries of weak acidic, neutral, and basic stationary phases for separating basic, neutral and acidic compounds of high polarity. Utilizing highest purity and enhanced mechanical stability silica and pure bonding reagents, Sepax HILIC bonded phases have been innovatively and specially designed to ensure maximum surface coverage, resulting in high stability of the stationary phases. The chemistry of monolayer formation is completely controlled that results in very reliable lot-to-lot and column-to-column reproducibility. The uniform, spherical Sepax HILIC particles have a nominal surface area of 300 m²/g with a controlled pore size of 120 Å. Sepax HILIC columns are packed with a proprietary slurry technique to achieve uniform and stable packing bed density for maximum column efficiency.

Structure of HILIC Stationary Phases

With the chemical structures shown in Figure 1, Polar-100 and Polar-Diol are neutral, polar phases, while Polar-Diol is more polar than Polar-100. Polar-Silica is a weak acidic phase. Polar-Pyridine and Polar-Imidazole are basic phases, while Imidazole phase is more basic than Pyridine phase.


Figure 1. Chemical structures of HILIC phases.

Characteristics

- HILIC phases with unique chemistries of acidic, neutral, and basic surfaces
- Ultra-pure silica particles with controlled pore sizes
- Available in particle size of 1.8, 2.2, 3, 5, and 10 μm
- High chemical stability for low leaching
- Available columns with ID in the range of 75 μm to 30 mm and length from 1 to 30 cm
- Available packings from grams to multi-Kilogram
- pH stability: 1.5 - 8.0
- Suitable for separations of polar pharmaceuticals, peptides, amino acids, and other compounds

Separation of Nucleosides


Figure 2. Separation of nucleosides by various HILIC columns (5μm, 4.6x150mm). Mobile phase: 90/10(V/V) = ACN/10mM ammonium acetate in H₂O. Flow rate: 1.0 mL/min. Ambient Temperature. UV 254 nm. Sample: 1. Uracil, 2. Adenosine, 3. Uridine, 4. Cytidine, 5. Guanosine.

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Sepax Technologies, Inc.